

FLAXSEED PILLOW FOR MICROWAVE USE

Instructions for use of the heat pillow

For your safety always follow these instructions.
Keep these instructions near your cushion.
Do not exceed recommended heating time and do not reheat until the cushion has returned to room temperature.
Allow to stand for 1 minute before applying to the body.

Heating:

Always read the instructions before use.
Always follow the instructions of the manufacturer of your microwave.
Do not use the cushion in a conventional oven or convector oven.
Carefully place the cushion in the microwave, making sure there is nothing inhibiting the movement of the turntable. If your microwave does not have a turntable, you should turn the cushion over halfway through the heating time to assure even heating. Always remove the shelf from the microwave before you warm up the cushion.
Take care when removing from microwave to avoid burns.

The cushion should feel warm not hot to touch.
Extreme overheating can lead to scorching or fire. Overheating can damage the cushion and its surroundings.

The microwave functions so that the heat slowly spreads from the inside to the outside. This means that it takes a little while for the cushion to feel hot on removal from the microwave. Do not leave in the microwave too long to avoid overheating!
Never leave the cushion unsupervised during heating.
Do not allow it to dry out.

Caution:

- Avoid contact with sensitive skin.
- Do not use on sleeping people.
- This item should not be used on babies or anyone unable to remove an inadvertently overheated pillow.
- Check the temperature of the pillow before using it.
- Do not use when symptoms of inflammation or infection are present.
- Persons with diabetes, circulatory, sensory, neurological and/or thermal problems should not use this product.
- Do not use in place of medical treatment.
- Store in a cool, dry place.
- To clean, wipe the sleeve with a damp cloth.
- Never wash the filling.
- Heating time 30 seconds intervals until desired warmth is achieved
- Maximum heating time should not exceed 90 seconds
- Make sure you consistently follow the times listed above and always check temperature of the pillow before using.
- Nature's choice is not responsible for any damage resulting from improper use or failure to follow the cushion's user instructions.
- Keep these instructions near your pillow.

Nature's choice b.v.
P.O. Box 626, 4870 AP Etten-Leur, The Netherlands
www.treets.com



ALWAYS READ THE INSTRUCTIONS BEFORE USE!

Fabric: 100% cotton
Filling: flaxseed