

PRECAUTIONS:

After the micturition (urination), it is recommended to change where URICLAK® is pressing on the penis, by either moving it ahead or backward.

It is recommended to use lubricant or moisturising milk while using URICLAK® to avoid possible irritation and pain in the skin.

You should replace your URICLAK® device every 6 months.

Do not use your URICLAK® device if the soft coating shows any form of deterioration.

Do not use the URICLAK® if the inner metal structure breaks.

URICLAK® should not be placed more than three hours pressing in the same place.

Not to use URICLAK® for more than eight hours a day.



URICLAK®

Distributed by:

DPM PROMOCIONES MÉDICAS, S.L.
Mare de Déu de la Salut, 26
08024 Barcelona (Spain)
Tel. +34 93 284 24 40
dpm@dpm-promedicas.com

**External device for male
urinary incontinence**

URICLAK® is an easy-to-use, transformable external controller for male urinary incontinence.

When wearing URICLAK®, use fingers to press either side to make it a circular ring that allows urination without having to remove the device.

It is practical, comfortable and easy to use.

This device avoids splashes or leaks onto clothes during urination as it does not need to be taken off or loosened like other devices. URICLAK® is very discreet and manageable.

It is very easy to change position of URICLAK® along the penis to avoid possible inconveniences.

URICLAK® is covered entirely with soft medical material so it is comfortable against the skin.

The innovative shape and inner moulding of the URICLAK® when closed, makes it better for blood flow while in use.

Using URICLAK® saves money on the use of absorbent compresses and pads.

INSTRUCTIONS:

With your thumb and forefinger, press on the narrow ends of your URICLAK®. Inside, the elongated bulge should be below and the two small bulges above (image 1).




Press URICLAK® inward (as per the pictures) making the device become circular. Place around the penis, centering it between the two outer bumps. Wear the device where it is most comfortable.



Do not press URICLAK® beyond the circular form so the metallic cross-bows double. This will cause them to lose their flexibility and pressure.

URICLAK® is supplied with four pairs of different sized pipes to use as regulators of pressure. Cut the desired pair from the plastic and insert into the circular holes at either side of the URICLAK®.

Choose which regulator pipe most suits your comfort and needs. Always start with the smallest diameter. A larger diameter will lower the pressure.



Important notice

To prolong the life of the URICLAK®, clean with a cotton swab slightly moistened with alcohol. Avoid contact with water.