

HairMax® Laser Device Model Comparisons



HAIRMAX®


| Features and Benefits | Laser Band 82 | Ultima 12 | Prima 7 |
|---------------------------------------|---------------|-----------|---------|
| Medical Grade Lasers | 82 | 12 | 7 |
| Approximate Treatment Time | 90 sec | 8 min | 15 min |
| Clinical Studies Proving Efficacy | √ | √ | √ |
| FDA Cleared for Males & Females* | √ | √ | √ |
| No Harmful Side Effects Ever Reported | √ | √ | √ |
| Cordless | √ | √ | √ |
| Patented Hair Parting Teeth | √ | √ | √ |
| Rechargeable Batteries | √ | √ | √ |
| 2-Year Warranty | √ | √ | √ |
| Worldwide Power Adapter | √ | √ | √ |

*HairMax laser devices are indicated to treat Androgenetic Alopecia, and promote hair growth in males who have Norwood Hamilton Classifications of IIa to V and in females who have Ludwig (Savin) I-4, II-1, II-2, or frontal patterns of hair loss and who both have Fitzpatrick Skin Types I to IV.

Are You a Candidate for HairMax® Laser Devices?


HairMax laser devices can be used by men and women with thinning hair or pattern baldness caused by a hereditary condition. Doctors use a system known as the Norwood Hamilton Classification (men) and the Ludwig-Savin Classification (women) to describe degrees of hair loss. Below are pictures representing different stages of hair loss. The shaded areas show the type of hair loss that can be treated effectively with the HairMax lasers devices. HairMax laser devices have not been clinically tested for hair loss outside of the shaded area.

Ludwig-Savin Classification - Females



HairMax Laser Devices have not been clinically tested for the hair loss outside of the shaded areas.

Norwood Hamilton Classification



HairMax laser devices were studied in light to brown skin tones (Types 1,2,3,4) as described in the Fitzpatrick Scale below. If your hair loss and skin tone fall into these categories, HairMax should work for you.

Fitzpatrick Skin Type Classification Guide: A way to describe your skin type based on your complexion and what happens to your skin in sunlight.

| | | |
|--|--------|---|
| Clinically Tested for these skin types | Type 1 | Very white or freckled skin, always burns with sun exposure (very fair; often in people with red or blond hair and blue eyes) |
| | Type 2 | White skin, usually burns with sun exposure (fair; often in people with red or blond hair and blue, green, or hazel eyes) |
| | Type 3 | White or olive skin tone, sometime burns with sun exposure (fair; seen in people with any hair or eye color) |
| | Type 4 | Brown skin, rarely burns with sun exposure (common in people of Mediterranean descent) |
| Not Clinically Tested for these skin types | Type 5 | Dark brown skin, very rarely burns with sun exposure (common in people of Middle-Eastern descent) |
| | Type 6 | Black skin, never burns with sun exposure |

*Adapted from Fitzpatrick's Dermatology in General Medicine. McGraw-Hill Professional; 5th edition.

Please note that HairMax was only tested on light to medium skin tones. This does not indicate that it will not work on darker skin tones.