

Tone and tighten your facial muscles and achieve a more youthful, contoured appearance by enhancing the conductivity of your Collagenius® EMS Facial Toner.

**Directions for use:**

Apply a small amount to the metal electrodes on your Collagenius® EMS Facial Toner before use. Wipe off with a clean damp cloth after use.

**First timers:**

Before using for the first time, apply a small amount to your forearm. If an allergic reaction occurs, discontinue use.

**Best before use/Batch number:**  
see package.

COLLAGENIUS®

TONING  
Conductive Gel

Increase conductivity,  
enhance your results

Manufactured for:  
Savantini Ltd.,  
Savantini House, Foster Street,  
Kingston-Upon-Hull,  
United Kingdom.  
HU8 8BT

**Ingredients:** Aqua, Glycerin,  
Hydroxyethylcellulose, Citric  
Acid, Potassium Sorbate,  
Sodium Benzoate

Keep out of reach of children.  
Store in a cool dark place.



  
see package



  
3261053



L13407\_1019

100ml