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---

## Light therapy unit BeautyLite BL 100

Successful use by acne and for Anti-Aging-Application by wrinkles, fine lines and blemishes

---

## Table of Contents-

Summary .....	9
Before operation .....	9
Safety instructions / Precautions .....	9
Contraindications .....	10
Explanation, applications, and benefits of the BeautyLite Therapy.....	10
What is the BeautyLite Therapy? .....	10
Applications of the BeautyLite Therapy.....	10
Benefits of using the BeautyLite BL 100.....	10
Operation of BeautyLite BL 100 .....	11
Basic adjustments .....	11
Different colours of LEDs for different applications .....	12
Cleaning, maintenance and disposal .....	12
Guarantee .....	12
Technical specifications and scope of delivery.....	12
Questions and answers .....	13

## Before operation

Dear customer,

We enjoy that you bought a product from DAVITA® and congratulate you to your choice. In this user manual particularly the functions and applications of the BeautyLite BL 100 are described in detail.

## Safety instructions / Precautions

Before using the unit the first time the user manual has to be read carefully. To reduce the risk of fire, burns, electric shock or other injury, before using your BeautyLite unit a careful study of the following important safety precautions is strongly recommended:

- ⚠ Warning – This Appliance must not be taken into the bathroom.
- ⚠ Do not immerse in water or any other liquid.
- ⚠ Do not use with wet hands.

Such situations above constitute a danger to life through electric shock.

- ⚠ Do not wrap the mains cord around the unit when storing.
- ⚠ This appliance should never be stored while still connected to the mains, always ensure the appliance is switched off, unplugged from the mains and that it has cooled down before storing.
- ⚠ The unit should not be left unattended when plugged in.
- ⚠ The unit should only be used for its intended use as described in this user manual.
- ⚠ Never continue to use this appliance if it is not working properly, or has been dropped, or if it has been dropped in water or any other liquid, or if it has a damaged plug or lead.

- ✚ Keep the mains lead away from heated surfaces.
- ✚ Do not stare into the light for too long.
- ✚ Do not repair the appliance yourself always consult an authorized technician.

## Contraindications

The BeautyLite should not be used by the following indications:

- ✚ Cancer
- ✚ Over any areas of infections
- ✚ Eczema
- ✚ Pregnant
- ✚ Epilepsy
- ✚ Cuts and open wounds
- ✚ Recent sun burn

Do not use around the eyes or mouth, or any sensitive parts of the face or body.

## Explanation, applications and benefits of the BeautyLite Therapy

### What is BeautyLite Therapy?

The BeautyLite Therapy is a hand held device which emits coloured light and is held over specific points located on the body. The light penetrates the skin to help improve the appearance of the skin and wellbeing. It is non-invasive, safe, and painless.

### Applications of the BeautyLite Therapy

BeautyLite Therapy is an innovative treatment method that uses pulses of red, yellow and blue light to stimulate and accelerate the body's natural healing process. The light stimulates the fibroblast cell activity in the skin, and results in an increase in the production of collagen and elastin. The skin will appear smoother, tighter and more radiant.

BeautyLite Therapy is a non-invasive treatment and is safe to use. BeautyLite Therapy does not heat the tissues the same way lasers do, the reason being that LED use longer wavelengths and therefore penetrate the tissues deeper, where lasers are more pinpointed in their delivery. The advantage of using LED is that you can treat a larger surface area in each treatment.

### Benefits of using the BeautyLite BL 100

The using of BeautyLite can:

- ✚ Significant reduce acne, wrinkles and fine lines
- ✚ Rejuvenating the skin
- ✚ Eliminate blemishes
- ✚ Improve irregular pigmentation
- ✚ Rejuvenate the skin
- ✚ Reduce redness and age spots
- ✚ Improve UV-damaged skin
- ✚ Increase the skin humidity
- ✚ Reduce the pore size
- ✚ Improve the appearance of skin
- ✚ Reduce the signs of aging

The results of the above help to improve the confidence as the skin will look better.

## Operation of BeautyLite BL 100

## Basic adjustments

- ✚ Before using ensure that the skin is completely clean and dry.
- ✚ Plug the mains adaptor into a mains socket, and the small Jack Plug (on the adaptor) fits into the back of the BeautyLite.
- ✚ Open the BeautyLite BL 100 by pressing the white button on the front of the unit down and the cover will be released. Lift the cover to open the unit.
- ✚ To switch "On" the BeautyLite unit, press the "ON" button.
- ✚ Once the „On“ button has been pressed the unit will „Beep“.
- ✚ On the LCD display the three primary colours and power levels are selected by default.
- ✚ The primary colours are Blue, Yellow and Red, and the secondary colours are Lilac (red & blue), Green (blue & yellow), Orange (red & yellow) Combined colours (blue, yellow & red).
- ✚ Select which colour you want to use by pressing the colour button or the combination as above. The unit will "beep" every time you press any of the buttons.
- ✚ Once you have selected your colour, select the power level by pressing the power level increase / decrease button.
- ✚ To begin your treatment push the Start button, this will „Beep“ and the lights will appear on the hand held unit.
- ✚ Gently run the lens of the BeautyLite unit slightly above the surface of the skin. Alternatively you can attach the enclosed plastic spacer for a more accurate reading. It is not essential to use the spacer.
- ✚ The BeautyLite unit programs are automatically set for 15 minutes.
- ✚ Treatment may be paused at anytime during the treatment by pressing the Pause button and the lights on the hand held unit will be switched off and the clock on the LCD will flash. To restart the treatment again press the Start button, the lights will come on again.
- ✚ If you want to stop the treatment at anytime press the Stop button, the lights will switch off and the display will reset.
- ✚ The BeautyLite BL 100 has an automatic „cut off“ switch after 5 minutes and the unit will have to be turned On again.
- ✚ When you have completed your treatment, please switch off the unit by pressing the „Off“ button and unplug the unit from the mains.
- ✚ After the BeautyLite unit has been put away safely, apply a small amount of the Soothing Moisturising lotion to the skin.

## Different colours of LEDs for different applications

### Blue light

Blue light has a calming effect on the body and also influences the elastic tissue content of new cell growth. But also targets the bacteria that causes acne and enhances collagen rejuvenation.

### Yellow light

Yellow light has a purifying effect on the toxin processing function whilst strengthening the endocrine functions. Firms and helps smooth the skin's surface and lighten scars. Helps with acne and wrinkles.

### Red light

Red light stimulates deep into the dermal layers and helps to activate the collagen and elastin production. Stimulates the natural process to rejuvenate skin, improves skin tone and texture, smoothing of fine lines and wrinkles.

### Lilac light (Red and blue light)

Red and blue light have an antibacterial effect and promote circulation. It can be used successfully against acne vulgaris. It can assist the healing and restores the skin's natural balance.

### Green light (Blue and yellow light)

Green light has a sedation and balancing effect on the mind and body and helps reduce inflammation and swelling of tissue.

### Orange light (Red and yellow light)

Orange light increases the blood and lymph circulation removing waste matter and helps nourish the cells.

### Blue, yellow and red light

The combination of the Red, Yellow, Blue light waves can help with rejuvenating the skin, helping heal acne and blemishes. Combined with the Yellow light it is both stimulation with the cells and effective healing on the surface. The combined colours will help clear away toxins within the cells.

## Cleaning, maintenance and disposal

A special maintenance is not necessary for this device. The unit should be switched off and unplugged before cleaning. To clean the unit a damp towel can be used. For cleaning the unit should never be immersed in water. Keep this product in a dry, cool and safe place. Do not use any harsh detergents, abrasives, solvents or cleaners.

Packaging: DAVITA® is joined to the system „Grüner Punkt“. All parts of the packaging are recyclable materials and environmentally acceptable to dispose. Units which are signed with the crossed wheellie bin symbol are not allowed to throw into the domestic waste. There is the commitment to dispose such electro- and Electronic "Altgeräte" separately. The communal public authority which is responsible for the disposal gives information regarding the possibilities of a regulated disposal.

## Guarantee

DAVITA® guarantees for his instrument for a period of 24 months after date of sale. Wear parts, such as batteries are excluded from the guarantee. The guarantee is valid for material defects and manufacturing faults and includes the provision of spare parts as well as labour. Guarantee expires if the instrument has been damaged, or used improperly or unauthorized intervention was carried out. Keep the purchase receipt as proof and submit it together with your claim for guarantee. DAVITA® will only fulfil the guarantee upon receipt of your purchase proof. Wear parts are generally excluded from the guarantee.

## Technical specifications and scope of delivery

Output:	6 V 300 mA
Input Voltage:	100 v / 240 V~, 50 – 60 Hz
Power:	1,8 W – 3,7 W
Weight:	0,51 Kg

Size of the system in the Beauty Case: 22 x 17 x 5 cm

Size of the hand unit: 15 x 7 x 4 cm

CE-mark, Art.No.: 20200

**Scope of delivery:**

BeautyLite BL 100 in the BeautyCase with integrated control unit incl. LCD-Display, hand held unit and mirror as well as adapter and detailed user manual.

## Questions and answers

**How does the BeautyLite work?**

Light enters the body as photons (which is energy) and then is absorbed by the photoreceptors within the cells. The photons strike the damaged tissue including skin, blood, muscle and bone which creates a cellular response and jump-starts the healing process.

BeautyLite works with visible light of different wave length between 460 – 655 nm. Light with a wave length of 650 – 750 nm will be recognized from the human eye as red, light with a wave length of 575 – 585 nm as yellow and light with a wave length of 420 – 490 nm as blue light. Every of these light colours has different effects on the human physis and psyche.

**How long will it take before results can be seen?**

This depends on what the individual is using the BeautyLite for. Generally there is an immediate response due to the enhanced blood circulation. Which in turn delivers an increase in the oxygen supply to the cells. It is especially successful in helping treat soft tissue injuries, leg ulcers, burns and acne.

**Is the BeautyLite painful to use?**

Not at all. Infact the sense of well being during and after using the BeautyLite Therapy is one of the benefits. There is only warmth and sense of relaxation, during and after use.

**How often should the BeautyLite be used?**

This depends on the severity of the problem being treated. The programs are set for 15 minutes and it is recommended the treatment should be carried out once a day to start with. When there is a visible difference it should be done approximately twice a week.

**What Time of the day should the BeautyLite be used?**

We believe that the best time of day is during the evening so that the cellular effect can continue for as long as possible whilst you sleep.

**Can anyone use the BeautyLite?**

Children under the age of 12 years should not use the BeautyLite unless it is recommended by the physician and only then under adult supervision. This is because their bodies are still in the stage of growth and children should not require treatments of this type at that stage of their life. Please check also the contraindications.

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DAVITA® Medizinische Produkte GmbH & Co. KG • Postfach 2004 • D-47518 Kleve