## 5 EASY STEPS

to get you up and running with your THERABATH PARAFFIN BATH!

**BEFORE YOU BEGIN**, locate and record the **Lot Number found on the bottom of your unit** (shown at right).

You will need this to register your warranty.

- 1007 03/14 MOEN USA R Medical Electronics Maple (RE) 188 Rates 100-128VAC. 1
- **Unpack** the **unit**; place it on a flat, secure surface.
- Pour up to 9 bags of refill paraffin wax in to the tank
- Place the plastic grille on top of the unmelted paraffin.
- Rest the lid on top of the grille. (the paraffin will compact as it melts)
- Attach plug; plug in the unit. (the paraffin will fully melt in 6 8 hours)







## **NOTE:**



During proper operation, the light turns on and off regularly (as the unit is heating), and is mostly off. The light is not an indication that the wax is melted.





## therabath quick tips

## cleaning

Clean unit when you see sediment has built up on bottom of inner tank. Unplug unit, remove lid, & grille and allow paraffin to completely harden

Plug in unit for approximately 10 minutes to melt edges of hard paraffin "cake" around the inside of the tank. Lay out paper towels to set wax "cake" on.

Turn unit upside down (keeping close to table surface). Holding the long edges of the tank, squeeze the tank outwards away from the wax, shaking gently until the cake releases and slides on to the

table surface. Blot up remaining paraffin in tank with a soft paper towel and wipe unit dry. HINT: a hair dryer will melt hard paraffin drips, then blot them up with a paper towel.

Replace fresh paraffin in the tank, place the grille on top of paraffin, and replace lid.

Plug in the unit.

After paraffin has melted (6-8 hours), add additional wax as necessary to bring up to required level.

www.therabath.com