

Collagenius Laser 180 Hair Growth Helmet

Safety Instructions

Before use please ensure you have read the following instructions;

- Do not look directly into the laser light. In built safety feature should ensure that lasers only switch on when the helmet is on your head
- Do not use when the hair is wet
- Expose as much of the scalp as possible before use.
- The helmet can not be used as a safety device to protect your head.
- Keep the device out of the reach of children
- Do not put it onto a fire or any area of high heat- rechargeable battery can rupture or explode
- Keep power cord and charging socket away from water
- Do not damage, heavily bend, or twist the power cord
- Only use the power cord provided

Using the Collagenius Laser 180 Hair Growth Helmet **should not** be used if you;

- Are a pregnant or breast-feeding woman
- Have a history of skin/scalp cancer
- Suffer from skin conditions such as eczema and psoriasis
- Suffer from a condition associated to allergies to light
- Sufferers of light induced epilepsy or chronic migraines
- Have high blood pressure, fever, heart conditions

Operating and Storage Environment

- Use the charger and lead supplied
- LED Laser Lights: 180 lights
- Helmet dimensions & weight: 266 x 194 x 133mm- 540g
- Laser Wavelength: 650nm
- Energy 5mW
- Illumination time: 25 mins
- Working voltage: DC 5V 1000mA
- Ensure the power is off, place in a dry and ventilated area and keep away from children