










# OSALIS

Health & Wellbeing



## INSTRUCTIONS

- 1 Plug the power cord in the input port on the back of your lamp.
- 2 Switch the lamp on by either pressing the power button  on the remote control or the back of the device.
- 3 To set the timer, press the timer button  before selecting your desired time. (10, 20, 30 or 60 mins)
- 4 To adjust brightness, press  to decrease or  to increase. Alternatively, select low, medium, or high brightness.
- 5 Use the    buttons to select the colour temperature.
- 6 To switch off the lamp, simply press the power button and unplug the power cord before storing your lamp away.

# LIGHT THERAPY FOR WELLBEING



## WARNINGS AND PRECAUTIONS

- Do not use with a power supply voltage other than 100-240 VAC.
- Do not operate device near water, flammable or combustible vapours, such as aerosol spray products, or where oxygen is being administered.
- Avoid placing in areas that are exposed to direct sunlight or close to heat radiating products such as heaters.
- Child supervision is recommended. This lamp contains small parts which could pose a choking hazard.
- Do not shorten or cut the power cord.
- Use with caution if you experience light sensitivity or have a pre-existing eye condition. Discontinue use if you experience discomfort. Certain medications (e.g. anti-depressants) may induce light sensitivity which can result in discomfort or hyperactivity. Consult with your healthcare provider before use.

---

**OSALIS**  
Health & Wellbeing



Savantini Ltd.,  
Savantini House, Foster Street,  
Kingston-Upon-Hull, HU8 8BT  
+44(0) 1482 496931  
[www.stressnomore.co.uk](http://www.stressnomore.co.uk)

IFU13821\_0322