

# Plan or Prevent Pregnancy Naturally?

## Compare our products...



### Time Level (Level 1)

Time, date, year	✓	✓	✓
Time zone and shift work	✓	✓	✓
Alarm time – Alarm clock	✓	✓	✓

### Fertility Level (Level 2)

Fertility display	✓	✓	✓
Display of ovulation day	✓	✓	✓
Data reading alternating with cycle day	✓	✓	✓
Data readings of previous days (in days)	250	180	99
Predicted future data readings (in days)	6	6	6

### Cycle Statistics (Lady-Comp® Level 3 – Lady-Comp Baby® Level 4)

Body Mass Index	✓	✓	✗
Number of stored cycles	✓	✓	✗
Average cycle length	✓	✓	✗
Average temperature increase after ovulation	✓	✓	✗
Average fluctuation of ovulation day	✓	✓	✗

### Planning Statistics (Lady-Comp Baby® Level 4)

Average length of high level	✓	✗	✗
Display of monophasic cycles in percent %	✓	✗	✗
Planning hits of the last 20 cycles	✓	✗	✗

### Pregnancy Test

Display of probable pregnancy after 15 days	✓	✓	✗
Display of confirmed pregnancy after 18 days	✓	✓	✗

### Other Functions

Computer self test and trial reading	✓	✓	✓
Measuring time window +/- 3 hours of preset alarm	✓	✓	✓
Menstruation input, automatic menstruation inquiry	✓	✓	✓
Input of past menstruation data (in days)	250	180	99
Cycle printout (in days)	250	180	99 - 180
Possible data deletion and input	✓	✓	✓
Possible upgrade with the planning statistics		✓	✗
Clinically tested Pearl Index	0.7	0.7	unter 1.0*
Product Life Expectancy (in years)	ca. 10	ca. 10	ca. 8-10
Extra costs and expendable items	Akku	Akku	Batterie
Size (in cm)	ø 14	ø 14	7 x 10
Weight (in g)	200	200	52

\*manufacturer information