



For additional information and your online order just visit our website www.saneo-line.com and www.facebook.de/saneoline



Massage (TENS/EMS)

by gently stimulating electrical impulses



effective, easy & safe as

- Relaxation massage
- Activation massage
- Flow massage
- Facial massage
- Leg-bum-tum massage
- Vein massage
- Circulation massage
- Tap massage
- Knead and grip massage
- Side of hand massage
- Vibration massage
- Hydro jet massage etc.



Further products SaneoLINE



I'm healthy!



Pain Relief

- chronic and acute pain
- headaches / migraines
- neck pains / tension
- pains in back and shoulder
- pains in joints
- sciatic nerve complaints
- rheumatic complaints and others



I'm in shape!



Muscle Training

- Firming of abs and buttocks
- Muscle growth/strengthening around the chest, arms and legs
- Pain relief for athletic injuries
- Warming up before training
- Relaxation after training
- Supporting active forms of sport

I'm relaxed!



tic Medizintechnik GmbH & Co. KG
 Endelner Feld 9 • D - 46286 Dorsten
 info@saneo-line.com
 Hotline +49 (0) 18 05 . 18 15 18
 (0,14 €/min. from a German landline, other providers may differ)

10/2013 / R3 CG

German Quality Brand



The **SaneoVITAL** TENS EMS massage machine allows you to enjoy a beneficial time-out after a strenuous or stressful day.

Depending on the massage program, **SaneoVITAL** can help you to relieve painful muscular tension, train your facial muscles or those in the stomach, legs and bottom regions. It can help to reduce pain caused by swelling, for example in the legs, or simply to relax perfectly, all using pleasant electrical impulses.

TENS EMS electrical stimulation has already been successfully employed by the medical profession for many years, due to its effectiveness having been proved by a number of scientific studies. Enjoy a heavenly massage of the upper or lower back, of the arms or legs - every part of the body needs revival occasionally.



With the SaneoVITAL TENS EMS massager, you can indulge yourself once in a while. Treat yourself!

- 17 massage programs
- 8 self-adhesive, reusable electrodes (50 x 50 mm, 32 mm round, 4 pieces each)
- 2 separately adjustable channels for treating several body regions



- easy-to-use
- incl. clip for fixing at belt
- large display for good readability
- timer function (automatic turning-off)
- comprehensive safety features
- high quality accessories (incl. 9V block battery)
- certified medical product



Technical characteristics

Digital 2 channel system
Output: 0 - 60 mA / channel, rectangular, charge compensated, biphasic pulses
Frequency: 2 - 100 Hz
Pulse width: 200 - 700 μ s

