

# Maxhim®

Intimate Essentials for Men

## **WATERPROOF VIBRATING PROSTATE MASSAGER**

.....  
Perfect shape for  
massaging the  
prostate

.....  
100% Waterproof

.....  
Powerful  
multi-speed  
vibration settings

## DIRECTIONS FOR USE

Insert 2 AAA batteries (not included) and wash in warm soapy water before use, making sure to rinse all soap off. Apply a condom, if preferred, for additional cleanliness. Followed by a water-based lubricant, such as the Maxhim Max Pro Therapeutic Prostate Massage Gel, to ease insertion. Insert the Maxhim Waterproof Vibrating Prostate Massager and adjust the vibration to a comfortable level, before beginning to gently stimulate the prostate. When finished, wash again before removing the batteries to extend their life.



## MAX PRO THERAPEUTIC MASSAGE GEL

Formulated especially for prostate massage

Water and silicone formula that won't dry out, ensuring maximum comfort  
Fragrance, colour and alcohol free

Buy from [www.stressnomore.co.uk](http://www.stressnomore.co.uk)

**Maxhim**<sup>®</sup>  
Intimate Essentials for Men

Savantini Ltd, Savantini  
House, Foster Street,  
Hull, HU8 8BT

---

This massager is for adult use as a non-medical preventative device. In the event of any medical condition or pain, discontinue use and consult a healthcare professional.

Made in China.

---

